



CellSonic

CellSonic Ltd of Switzerland

Registered number CH-217-3540042-9
 Focused and Radial Shockwave Lithotripters
www.cellsonic-medical.com



CellSonic (electro-hydraulic shockwaves)
 Focused shockwaves from 5 mm to infinity
 Each head 50,000 shocks
 Supplied with four heads
 Variable power: wounds, bones and sports injuries
 Ten year warranty

Typical cost of heel spur treatment:	
Machine amortised over 10 years	\$2
1000 shocks	\$7
Total	\$9
Typical cost of wound treatment:	
Machine amortised over 10 years	\$2
3000 shocks	\$21
Total	\$23



SynchroWave radial shockwave machine for physiotherapists. With electro-therapy stimulation synchronised to the shockwave. Frequencies and powers adjustable. External compressor required.



CellSonic Gel
for perfect transmission of shockwaves.
With aloe vera for skin care.



CellSonic Physio
Radial shockwaves for physiotherapy
Adjustable frequency and power
Maintenance free compressor
Amortise over 4 years for machine cost of \$1-50 hour.

LITHOSPLIT 4.1 Digital
ISWT lithotripter for safe, sure and cost effective kidney stone removal.
Runs on hospital oxygen supply.
Supplied in carrying case with all you need to operate; probes, connecting tubes, foot switch and tools.
Range of accessories include bridge for bladder stones, adapter for ureteroscope and suction assembly for PCNL.



Lithosplit high precision Handpiece
accepts various probe sizes.

Manufactured in a factory approved to ISO 9001:2001 & ISO 13485:2003 for the design and manufacture of medical equipments

We also make
Uroflometry Systems
Insufflators
Cold Light Fountains
Various dermatology machines
and many more

Please check our website for full details.
www.cellsonic-medical.com
Your local distributor will help you
If you do not find a distributor in your area,
please contact us at
sales@cellsonic-medical.com